

SUPPLEMENTARY INFORMATION

1

2

MATERIALS AND METHODS

MTS assay

3 The cytotoxicity of sigma-2 ligands was determined using the CellTiter96 Aqueous One
4 Solution Assay (Promega, Madison, WI). Briefly, EMT-6 or MDA-MB-435 cells were
5 seeded in a 96-well plate at a density of 2000 cells/well on the day prior to treatment with
6 sigma-2 selective ligands. After a 24 or 48 hour treatment, the CellTiter 96 AQueous One
7 Solution Reagent was added to each well, and the plate incubated for 2 hours at 37°C.
8 The plate was then read at 490 nm in a Victor³ plate reader (PerkinElmer Life and
9 Analytical Sciences, Shelton, CT). The EC₅₀ value, defined as the concentration of the
10 sigma ligand required to inhibit cell viability by 50% relative to untreated cells, was
11 determined from the dose response curve for each cell line.
12
13

14

LDH assay

15 The lactate dehydrogenase (LDH) release assay was performed using the Cytotox 96®
16 Non-Radioactive Cytotoxicity Assay (Promega, Madison, WI) according to the
17 manufacturer's protocol. Briefly, MDA-MB-435 cells were seeded at 2×10^3 cells/well in
18 96-well plates on the day prior to treatment with the sigma-2 ligands. After a 24 hour
19 treatment with the various sigma ligands, plates were centrifuged at $250 \times g$ for 5
20 minutes. 50 µl of supernatant was removed and transferred to a new plate. 50 µl of the
21 reconstituted substrate mix was then added to each well and the plate was allowed to
22

1 incubate at room temperature for 30 minutes, protected from light. After 30 minutes, 50
2 µl of stop solution was added to each well and the plate was read at 490 nm in a Victor³
3 plate reader (PerkinElmer Life and Analytical Sciences, Shelton, CT).

4

5 **Detection of intracellular caspase-3 activity**

6 The activation of endogenous caspase-3 by the sigma ligands was measured using the
7 CellProbe HT caspase-3 whole cell assay (Beckman Coulter, Fullerton, CA). Briefly, the
8 cell permeable, nonfluorescent bisamide substrate, Z-DEVD-R110, is added to the intact
9 cells to detect apoptosis. Upon specific cleavage by caspase-3, Z-DEVD-R110 is
10 converted to a fluorescent compound whose fluorescent signal is proportional to the
11 amount of caspase-3 activity in the cells. After adding the substrate Z-DEVD-R110, the
12 plate is incubated for 1.0-1.5 hour, and the resulting fluorescence measured using a
13 Victor³ microplate fluorometer (PerkinElmer Life and Analytical Sciences, Shelton, CT)
14 at excitation and emission wavelengths of 485 nm and 535 nm, respectively.

15

16 **The annexin V assay**

17 Annexin V positive cells were quantified with an annexin V-FITC assay (Apoptosis
18 Detection Kit, R&D Systems, Minneapolis, MN) and quantified using flow cytometry.
19 Annexin V binds to phosphatidyl serine, which is translocated from the inner leaflet of
20 the plasma membrane to the outer leaflet in the early stage of apoptosis.

21

22

1 **SUPPLEMENTARY FIGURE LENGENDS**

2 **Supplementary Figure 1** Sigma-2 ligands increased the percentage of annexin V and
3 TUNEL positive cells. EMT-6 cells were either left untreated or treated with WC-26 (40
4 μM) for 48 h, SV119 (100 μM) for 16 h, or RHM-138 (40 μM) for 16 h; annexin V
5 positive cells and TUNEL positive cells were quantified by flow cytometry. The data
6 showed that the percentage of annexin V positive cells was increased from $7 \pm 1.1\%$ for
7 the untreated cells to $96.6 \pm 2.7\%$, $73.5 \pm 2.4\%$, and $48 \pm 6.4\%$ after treatment with WC-
8 26, SV119 and RHM-138, respectively. TUNEL positive cells increased from $4.57 \pm$
9 2.4% for the untreated cells to $48 \pm 2.4\%$, $36 \pm 3.0\%$, and $39.2 \pm 3.5\%$ for the cells
10 treated with WC-26, SV119 and RHM-138, respectively. $*p < 0.0005$ compared to the
11 untreated control.

12

13 **Supplementary Figure 2** RHM-138 induced caspase 8 and 9 activation. MDA-MB435
14 cells were treated with 40 μM RHM-138 for 0, 1, 2, 4, 8, 16 and 24 h. Caspase 8 and 9
15 cleavage was analyzed by western blot. RHM-138 induced caspase 8 and 9 cleavage,
16 suggesting that RHM-138 may trigger apoptosis via both the intrinsic and extrinsic
17 pathways.

18