

S2 Table. Tests for non-specific activity of compounds 1, 2 and 3

	% change in rate		
	compound 1	compound 2	compound 3
10x enzyme*			
1 nM TEM	+65%	+38%	-47%
10 nM TEM	+31%	+33%	-48%
Detergent-dependence*			
0% Triton-X	+2%	+30%	-40%
0.01% Triton-X	+66%	+33%	-45%
Alternative enzymes			
chymotrypsin	-5%	-1%	-2%
β-galactosidase	+2%	+5%	+3%

*Rates measured at 50 μ M nitrocefin.